

**Table 1. Multifocal Soft Contact Lenses on the Market\*\***

Multifocal Lens Name	Manufacturer	Material	Diameter (mm) / H <sub>2</sub> O%	BC	DK / CT / OZ (mm)
Concise*	ABB Optical	Polymacon C	14.5 / 38%	8.4, 8.6	
Concise MTO*	ABB Optical	Polymacon C	13.0 to 16.0 in 0.10 steps / 38%	7.8 to 10.0 in 0.10 steps	
Definitive MTO*	ABB Optical	Efrolfilcon A	13.5 to 15.5 in 0.10 steps / 74%	7.5 to 9.6 in 0.30 steps	
Definitive*	ABB Optical	Efrolfilcon A	14.2, 14.5 / 74%	8.3, 8.6, 8.9	
Ocu-Flex-53*	ABB Optical	Ocuflcon B	14.5 / 53%	8.4, 8.6	18.0 / varies / 1.9, 2.18
Air Optix Aqua	Alcon	Lotrafilcon B	14.2 / 33%	8.6	110, 138 @ -3.00D / 0.08 mm @ -3.00D / 0.10 @ -3.00D
Focus Dailies Progressives	Alcon	Nelfilcon A	13.8 / 69%	8.6	0.10 @ -3.00D
Dailies AquaComfort Plus	Alcon	Nelfilcon A	14.0 / 69%	8.7	26.0 / 0.10 (-3.00) / varies
Dailies Total1 Multifocal	Alcon	Delefilcon A	14.1 / 33%	8.5	156 @ -3.00D / 0.09
Astera	Alden Optical	Hioxifilcon D	10.0 to 16.0 (14.5 std.) / 54%	6.5 to 9.7 (8.3, 8.6, 8.9 std.)	21 / 0.1 / 8.0
Intelliwave Quarterly Replacement*	Art Optical Contact Lens	Efrolfilcon A	74%	MTO	60 / varies / 8.4
Intelliwave Semi-Annual Replacement*	Art Optical Contact Lens	Acofilcon A, Hioxifilcon B	48%, 49%	MTO	16 / varies / 8.4
Biotrue OneDay for Presbyopia	Bausch + Lomb	Nesofilcon A	14.2 / 78%	8.6	42 @ -300D / 0.10mm @ -3.00D / 9.0mm @ -3.00 D
PureVision 2 for Presbyopia	Bausch + Lomb	Balafilcon A	14.0 / 36%	8.6	130 @ -3.00D / 0.07mm @ -3.00D / 9.00mm @ -3.00D
PureVision	Bausch + Lomb	Balafilcon A	14.0 / 36%	8.6	101 @ -3.00D / 0.09mm @ -3.00D / 8.0mm @ -3.00D
SofLens	Bausch + Lomb	Balafilcon A	14.5 / 38%	8.5, 8.8	8.4 @ -3.00D / 0.10mm @ -3.00D / 8.0mm @ -3.00D
Ultra for Presbyopia	Bausch + Lomb	Samfilcon A	14.2 / 46%	8.5	163 Dk/t @ center for -3.00D / 0.07mm for -3.00D
Essential Soft Toric	Blanchard Contact Lens	Hioxifilcon B	14.2, 14.4, 14.8 / 48%	14.2 to 8.2, 14.4 to 8.5, 14.8 to 8.8	15 / 0.12 / N/A
Esstech PS	Blanchard Contact Lens	Polymacon	14.5 / 38.6%	8.3, 8.7, 9.1	8.4 / 0.14 (-), 0.15 (+) / N/A
Esstech PSD	Blanchard Contact Lens	Polymacon	14.0 / 38%	8.3, 8.7, 9.1	8.4 / 0.14 (-), 0.15 (+) / N/A
Quattro	Blanchard Contact Lens	Hioxifilcon B	14.2 / 48%	14.2 to 8.4, 14.5 to 8.8	15 / 0.08 (-), 0.12 (+) / N/A
CO Soft 55 Crescent Custom Bi-focal*	California Optics	Methafilcon A	14.5, 15.0 / 55%	8.6, 8.9, 9.2	18.8 / varies
CO Soft 55 Custom Progressive*	California Optics	Methafilcon A	14.5, 15.0 / 55%	8.3, 8.6, 8.9, 9.2	18.8 / varies
Biofinity	CooperVision	Silicone Hydrogel	14	8.6	128 / 0.08 (-3.00)
Clarity 1-day	CooperVision	Somofilcon A	14.1 / 56%	8.6	60 / 0.07
Proclear 1-day	CooperVision	Omafilcon A	14.2 / 60%	8.7	25 / 0.09 (-3.00)
Proclear EP	CooperVision	Omafilcon A	14.4 / 60%	8.7	33 / 0.16 (-3.00)
Proclear*	CooperVision	Omafilcon A	14.4 / 62%	8.7	34 / 0.16
Proclear XR	CooperVision	Omafilcon A	14.4 / 59%	8.7	34 / 0.16 (-3.00)
Triton Soft Translating Bi-Focal	GelFlex USA	Hioxifilcon B	15.0, 13.4, 14.5 to 13.4 / 48%	8.0, 8.2, 8.4, 8.6, 8.8, 9.0, 9.2	Varies / 0.15 / varies
iSight MCL	GP Specialists	Hioxifilcon A	13.0 to 15.5	8.0 to 9.5	59.8 / 0.08 / 1.8 to 3.5
1-Day Acuvue Moist	Johnson & Johnson	Etafilcon A	14.3 / 43%	8.4	26.0 / 0.084 (-3.00)
Acuvue Oasys for Presbyopia	Johnson & Johnson	Senofilcon A	14.3 / 38%	8.4	103 / 0.070 (-3.00)
Metrofocal*	Metro Optics	Polymacon	14.0 / 38%	8.6, 8.9	8.4 / 0.10 / 9.0
Definition Aberration Control	Optical Connection	Methafilcon A	14.2 / 55%	8.6	19.5 / 0.09 @ -3.00 / 11.0 @ -3.00
54 Bifocal/Multifocal Sphere*	SpecialEyes	Hioxifilcon D	12.5 to 16.0 in 0.1mm steps / 54%	6.9 to 9.5 in 0.1mm steps	23 / 0.10 to 0.30 / Custom Distance and Near
Duette	SynergEyes	Perfilcon A	14.5	7.1 to 8.3 in 0.2mm steps	130
SynergEyes PS	SynergEyes	Pafufocon D hem-iberfilcon A	14.5 / 27%	7.20 to 8.2 in 0.1mm steps	100 / 0.23mm, -3.00D / 7.8
C Vue 55*	Unilens	Methafilcon A	14.4 / 55%	8.5, 8.8	18.8 / varies / varies
C Vue Advanced HydraVue*	Unilens	Efrolfilcon A	13.5 to 16.0 / 74%	8.0 to 9.5	60 / 0.12 to 3.00 / 8.0
C Vue Advanced*	Unilens	Hioxifilcon D	13.5 to 16.0 / 55%	8.0 to 9.5	23 / 0.07 + 0.40 / 8.0
C Vue	Unilens	Polymacon	14.5 / 36%	8.5, 8.8	8.4 / 0.10 / 8.0
C Vue ADDvantage	Unilens	Balafilcon A	14.0 / 36%	8.6	91 / 0.07 to 3.00 / 9.0
C Vue HydraVue	Unilens	Balafilcon A	14 / 36%	8.6	91 / 0.09 to 3.00 / 9.0
LifeStyle 4-Vue	Unilens	Polymacon	14.5 / 38%	8.8	16 / 0.17
LifeStyle 4-Vue Hi-Add	Unilens	Polymacon	14.5 / 38%	8.5, 8.8	16 / varies / varies
LifeStyle MV2*	Unilens	Polymacon	14.5 / 38%	8.5, 8.8	16 / varies
LifeStyle Xtra	Unilens	Polymacon	14.5 / 38%	8.5	16 / varies / varies
SimulVue	Unilens	Hefilcon A	14.5 / 45%	8.4, 8.7, 9.0	11.6 / 0.14 / 9.0
SimulVue 38	Unilens	Polymacon	14.5 / 45%	8.4, 8.7, 9.0	8.4 / 0.10 / 8.0
Softsite	Unilens	Polymacon	14.5 / 38%	8.4, 8.7	8.4 / 0.10 / 8.0
Unilens	Unilens	Hefilcon A	14.5 / 45%	8.4, 8.7, 9.0	11.6 / 0.16 / 9.0
Unilens 38	Unilens	Polymacon	14.5 / 38%	8.4, 8.7, 9.0	8.4 / 0.10 / 8.0
Unilens EMA	Unilens	Polymacon	14.5 / 38%	8.5, 8.8	8.4 / 0.10 / 8.0
UCL Bifocal*	United Contact Lens	Ocuflcon C	14.5, 15.0 / 55%	8.3, 8.6, 8.9, 9.2	18.8 / 0.08, 0.30 / 1.7, 2.0, 2.3, 2.6, 2.9, 3.2
UCL Multifocal*	United Contact Lens	Ocuflcon C	14.7 / 55%	8.3, 8.7, 9.0	18.8 / 0.06 to 0.22 / 8.0 to 10.0
Horizon Bicon Sphere*	X-Cel Specialty Contacts	Various	14.0 to 15.0, 1mm steps / 49, 55, 59, 74%	8.0 to 9.2 in 0.1mm steps	0.14 to 0.37, varies / 1.8, 2.0, 2.5, 3.0, 3.5, 4.0
Horizon Progressive*	X-Cel Specialty Contacts	Various	14.0 to 15.0 in 0.1mm steps / 49, 55, 59%	8.0 to 9.2 in 0.1mm steps	varies
IdealSoft PS	X-Cel Specialty Contacts	Polymacon	14.0 / 38%	8.3, 8.6, 8.9	8.4 / varies

\*Toric versions also available. \*\*Source: 2016 *Review of Cornea & Contact Lenses* Annual Lens Guide. Data is provided by manufacturers and is not independently verified.